
MODE D'EMPLOI

Electrode de REF. : OGR002 - AR02318

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
OPERATING INSTRUCTIONS

REF. Electrode : OGR002 - AR02318

The AR02318 electrode is a Mercury Oxide (Hg/HgO) 120 mm ROD reference electrode. It is ideal for measurements in alkaline environments.

This electrode contains mercury. Dispose according to Local, State or Federal Regulations.

T+: very toxic

N: Dangerous for the environment

Do not throw in the bin

Electrochem
OrigaLys

Preparation for measurement:

1 Remove the protection cap from the electrode and any seals covering the filling hole and the porous plug (pin).

2 Before starting a measurement, remove the clip which closes the electrode filling hole. Remember to replace the clip at the end of measurements.

3 Check the level of the filling solution. It should be approximately 0.5 cm below the filling hole.

For tests and on delivery, the electrodes are filled with 0.1 M potassium hydroxide.

However, we highly recommend that you replace this solution with a solution which has a better stability, e.g. KOH 1M.

If carbonates are formed in the KOH solution, the reference potential will drift. Other concentrations or other hydroxides may therefore be used.

Note: extreme care must be taken when working with strong concentrations of potassium hydroxide as they are very corrosive. In case of contact with the skin, wash immediately with plenty of water.

4 In order to remove air bubbles trapped inside the electrode, shake the electrode holding it at its head with the porous pin down.

Specifications:

- Temperature range: -5 à 50°C

Note: the reference potential at 25°C depends on the type of filling solution used:

- For an electrode filled with 1M KOH it is approximately -129 mV.
- For an electrode filled with 0.1M KOH it is approximately - 70 mV, with respect to a saturated calomel electrode.

N.B.: In case of contact with mercury, wash exposed areas immediately with plenty of water.

Maintenance:

1 Electrode contamination is a major cause of faulty measurements.

2 The electrode should be rinsed with distilled water after measurements.

3 Check frequently the level of filling solution.

In case of deposits which cover the electrode, clean the electrode with:

- 4
- a solution of acid (0.1M HCl, 0.1M HNO₃): mineral salt deposits, etc.
 - Normal Cleaning Solution : greasy or oily deposits.

Note : the porous surface of the electrode can be cleaned using very fine abrasive paper.

Storage:

- Between measurements: leave it in a potassium hydroxide solution.
- Several days: cover the porous pin or fibres with paraffin film or with the electrode clip and immerse the electrode in potassium hydroxide solution.
- Long-term storage: it is recommended to empty the electrode and refill it with a fresh solution of KOH 0.1M. Do not forget to fill the protection cap with paraffin film and place the electrode in its box and store at room temperature.

Accessories:

- AR01210: Electrode Cable, with BNC plug, 1m
- AR01211: Electrode Cable, with BNC plug, 0,16m
- AR01206: Electrode Cable, with Banana Ø 4mm plug, 1m
- Normal Cleaning solution 150 ml
- Normal Storage solution 150 ml
- Saturated KCL Solution 150 ml