MODE D'EMPLOI

Electrode de REF. Calomel : OGR004 - AR02320

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
OPERATING INSTRUCTIONS

REF. Calomel Electrode : OGR004 - AR02320

Electrode AR02320 is a 103 mm Saturated Calomel Reference Electrode (ECS).

This electrode contains mercury. Dispose according to Local, State or Federal Regulations.

T4 : Very toxic

N : Dangerous for the environment

Do not throw in the bin

Preparation for Measurement:

- 1 Remove the protection cap from the electrode and any seals covering the filling hole and the porous plug (pin).
- 2 Before starting a measurement, remove the clip which closes the electrode filling hole. Remember to replace the clip at the end of measurements
- 3 Check the level of the filling solution. It should be approximately 0.5 cm below the filling hole.
- 4 In order to remove air bubbles trapped inside the electrode, shake the electrode holding it at its head with the porous pin down.

Specifications:

- Temperature range: -5 to 50°C

N.B. : In case of contact with mercury, wash exposed areas immediately with plenty of water.

Accessories:

- AR01210 : Electrode Cable, with BNC plug, 1m
- AR01211 : Electrode Cable, with BNC plug, 0,16m
- AR01206 : Electrode Cable, with Banana Ø 4mm plug, 1m
- Normal Cleaning solution 150 ml
- Normal Storage solution 150 ml
- Saturated KCL Solution 150 ml

Maintenance:

- 1 Electrode contamination is a major cause of faulty measurements.
- 2 The electrode should be rinsed with distilled water after measurements.
- 3 Check frequently the level of filling solution.
- 4 In case of deposits which cover the electrode, clean the electrode with:
 - a solution of acid (0.1M HCl, 0.1M HNO₃) : mineral salt deposits, etc.
 - Normal Cleaning Solution : greasy or oily deposits.

Note : The porous plug of the electrode can be cleaned using a fine abrasive paper.

Storage:

- Between measurements: leave it in a Saturated KCl Solution.
- Several days: seal the filling hole with paraffin film or with the electrode clip and immerse the electrode in a Saturated KCl Solution using the accompanying cap.