
MODE D'EMPLOI

Electrode de REF. Ag/AgCl : OGR006 - AR02322

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
OPERATING INSTRUCTIONS

REF. Ag/AgCl Electrode: OGR006 - AR02322

The AR02322 electrode is a 120 mm silver chloride (Ag/AgCl) reference electrode.

electrochem
OrigaLys

Preparation for Measurement:

- 1 Remove the film covering the filling port and the cap protecting the porous tube.
- 2 Remove the clip obstructing the filling port before each series of measurements. This clip will be replaced at the end of the period of use.
- 3 Check the level of the reference electrode filling solution: it should be approximately 0.5 cm below the filling hole. Refill, if necessary, with a 3M KCl solution saturated with AgCl. The concentration of KCl can be varied if desired, but it is advisable to use a high concentration. Remember to saturate the solution with AgCl. For applications in acetic acid: Empty the electrode and fill it with a saturated solution of KCl in acetic acid. Saturate the solution with AgCl.
- 4 To eliminate any air bubbles that may have formed inside the electrode, shake the electrode like a thermometer.

Specifications:

- Temperature range: -5 à 50°C

Accessories:

- AR01210 : Electrode Cable, with BNC plug, 1m
- AR01211 : Electrode Cable, with BNC plug, 0.25m
- AR01206 : Electrode Cable, with Banana Ø 4mm plug, 1m
- Normal cleaning solution 150 ml
- Normal storage solution 150 ml
- Saturated KCl solution
- 3M KCl solution saturated with AgCl, 150 ml

Maintenance:

- 1 Electrode contamination is a major cause of faulty measurements.
- 2 The electrode should be rinsed with distilled water after measurements.
- 3 Check frequently the level of filling solution.
- 4 In case of deposits which cover the electrode, clean the electrode with:
 - A solution of acid (0.1M HCl, 0.1M HNO₃): mineral salt deposits, etc.
 - Normal Cleaning Solution : greasy or oily deposits

Note: The porous plug of the electrode can be cleaned using a fine abrasive paper

Storage:

- Between measurements: leave it in 3M KCl solution.
- Overnight or longer: seal the filling hole with paraffin film or with the electrode clip and fit back in place the protection cap filled with the 3M KCl Solution.